


Worm Club Journal Club

Joey Doll
August 1, 2011



Automated on-chip rapid microscopy, phenotyping and sorting of *C. elegans*





Kwanghun Chung^{1,3}, Matthew M Crane^{2,3} & Hang Lu^{1,2}

NATURE METHODS | VOL.5 NO.7 | JULY 2008 | 637

Noncontact microrheology at acoustic frequencies using frequency-modulated atomic force microscopy

Núria Gavara & Richard S Chadwick

650 | VOL.7 NO.8 | AUGUST 2010 | **NATURE METHODS**



FACS for worms

Sort based upon complex reporter expression patterns.

Two layer PDMS device

- worm layer
- valve/cooling layer


Important features...

- simple loading system
- reproducible worm placement
- worm cooling (vs. drugs/clamp)
- works on any microscope
- all features $>20\mu\text{m}$


900 worms/hour (epi)

150 worms/hour (confocal)


a Worm entering







b Worm being positioned




c Worm being cooled and imaged




d Worm exiting




Frame 1: red
Wait 270 msec
Frame 2: green



Mechanical clamp





4C cooling





Confocal z-stack to determine number of AWC neurons expressing *pstr-2-gpf*.



Focal plane 1



Focal plane 2




GFP intensity




Thresholded





k








Measure soft material mechanics (10-200 kPa) based upon hydrodynamic forces on an AFM cantilever.




Developed at NIDCD lab to measure tectorial and basilar membrane mechanics at acoustic frequencies.



- Measure resonant frequency (90 degree phase) as a function of sample-bead gap.
- The stiffness of the squeezed film depends on the sample mechanics, and affects freq.






Bead embedded in 6 kPa gel

Contact mode F-D curve
(Forces ~ 1000 pN)

Oscillating bead
(Forces ~ 1 pN)



Bead measurements perturb the sample less than the conventional F-D curves.




Table 1 | Viscoelastic properties of the tectorial membrane

Location	E (kPa)	μ_{eff} (Pa \times s)	G' (kPa)	G'' (kPa)	Loss tangent (G''/G')
Apex	3.3 ± 2.3	0.007 ± 0.006	1.1 ± 0.8	1.3 ± 1.1	1.2 ± 1.9
Middle	5.7 ± 3.3	0.01 ± 0.008	1.9 ± 1.1	1.8 ± 1.5	0.9 ± 1.3
Base	14 ± 16	0.03 ± 0.04	4.7 ± 5.3	5.6 ± 7.4	1.2 ± 2.9

Data are mean \pm s.e.m. ($n = 4$). One-way ANOVA reported significant differences in E ($P = 0.006$) and μ_{eff} ($P = 0.007$) along the length of the cochlea.

