

1 mm



Worm Club
Update





Joey Doll
3/29/10

The Plan

- Introduction
- Problems and solutions
 - Curvature induced stiffening
 - Accidental resistor etching
 - Actuator to sensor crosstalk
 - Experimental issues
- Current work
 - New die layout
 - Design specs
- What's next


The Devices

Piezoresistive
force sensor
(PR)




Piezoresistive force sensor
with integrated
piezoelectric actuation
(PR+PE)

Applications of a Faster Actuator and Sensor




Sense of touch (TRNs)



Single Molecule Force Spectroscopy


Images from Worm Atlas, <http://www.scienceinschool.org/print/210> and
Clues to the cochlear amplifier from the turtle ear, Fettiplace, Ricci and Hackney, 2001



Sense of hearing (hair cells)

Performance Goals



- Actuator
 - Rise time = 1-10 microseconds
 - Tip deflection < 1 micron
- Sensor
- Closed loop control
- Operation in fluid



Where Things Stand


- Made first generation PR and PR+PE devices
- Obtained lots of useful characterization data
- The fabrication process mostly works
- But some issues...
 - Force offsets
 - High resistances
 - Crosstalk
- Designing the next round of devices now

Curvature Induced Stiffening




- A large force is required to overcome the “potato chip” effect and start bending the PR
- Due to surface stresses (probably dopants) and the way the cantilever is designed

Solving the Stiffening Problem



- 1) Add a stress compensation layer
Blue = dopants (tensile)
Orange = surface oxide (compressive)





- 2) Modify the cantilever design to reduce stress sensitivity (eliminate long transverse sections)

- 3) Measure film stresses whenever possible and choose thicknesses accordingly

Dopant concentration stress model based upon experimental data and modeling




Accidental Resistor Etching



Problem: Resistances were about 10x higher than expected and depended heavily on voltage

Solution: Protect the piezoresistor during AlN deposition. Simple fab solution.

Experimental Crosstalk Data




Crosstalk is important due to the signal magnitude difference:

PR output = microvolts before amplification


PE drive = volts

Ideally want $< 1 \mu\text{V/V}$ crosstalk at 100 kHz (40,000x improvement from current)

Crosstalk Schematic



- Aluminum
- Silicon Dioxide
- Titanium
- Aluminum Nitride
- Silicon




Crosstalk Model

$$\frac{V_A}{V_{drive}} \approx \frac{Z_{Ti,Bottom}}{Z_{AlN,Top}}$$

$$\frac{V_B}{V_A} \approx \frac{Z_{LTO}}{2Z_{AlN,Bottom}}$$

$$\frac{V_B}{V_{drive}} \approx \frac{Z_{Ti,Bottom} Z_{LTO}}{2Z_{AlN,Bottom} Z_{AlN,Top}}$$





- Model predicts 3 mV/V @ 100 kHz (vs 40 mV/V measured)
- Close enough given the uncertainties involved

Crosstalk Solutions


- Substrate ground
 - Removes the LTO from the equation on the last slide → 30,000x improvement
 - Will probably use a p-n junction too for further improvement
- Symmetric device and circuit layout
 - Right now the crosstalk only goes into one input of the instrumentation amplifier
 - Utilize common mode rejection of amplifier to attenuate anything that gets around the isolation well
- >40,000x reduction seems reasonable
 - At least -80 dB from substrate ground
 - -100 dB from CMRR (with perfect device layout)

Some Experimental Issues

- Parylene coating
 - Thin (< 200 nm) coating needed to passivate but not affect mechanics, fixed
- Bead adhesion to parylene
 - Fixed, was just an issue of using the right epoxy + curing
- Wirebond passivation
 - The new, thin parylene doesn't cover the bondpad wires perfectly, will start encasing in epoxy
- Vibrations
 - ~ 100 nm vibrations on current setup
 - Moving to our AFM, want < 1 nm vibration
- Big beads on little cantilevers
 - Need to include bead mass in design optimization




New Device Layout



Old

1 mm



New

Substrate
contact

Temp. Comp.


Force Meas.

Thick extender


All resistors on die
(for crosstalk purposes)

New Device Layout

40 degrees



10 degrees



- This is for very, very long cantilevers (800 μm)
- The difference is much larger for short (high freq) levers

The Next Designs

- Combining the PE and PR optimization code is trickier than I thought – no designs ready yet
- General improvements
 - Shift cantilever to one side
 - Make it easy to use at an angle
 - Add die labeling
 - Define cantilever from the frontside
 - Skinnier PCB tip

Double-Checking the Specs

- Mammalian hair cells
 - 1 microsecond rise time → 300 kHz resonant frequency
 - Current = 10 microseconds (30 kHz)
 - 50 pN force resolution (plausible but hard @ 300 kHz)
 - Max deflection > 1 micron
 - Inner: 7-10 micron diameter bead
 - Outer: 3-5 micron diameter bead
- Turtle hair cells
 - Stiffness = 2-10 mN/m
 - 2-5 micron diameter bead
- General
 - Soft probes (mechanics) → ~1 mN/m
 - Stiff probes (kinetics) → much stiffer than the soft

What's Next

- Finish characterizing PR devices and finalize models
- Modify AFM for PR measurements and perform demo streptavidin-biotin experiments
- Extend optimization code to PR+PE designs
- Generate + review final designs w/ interested people
- Fab run



Thank You