




Research Update:


**Nonlinearities and
exploding devices**

Joey Doll
2/10/2010



Motivation



Sense of touch (TRNs)



Single Molecule Force Spectroscopy (SMFS)



Sense of hearing (hair cells)



Performance Goals

- Rise time
 - 10 microseconds (turtle hair cell)
 - 1 microsecond (mammalian hair cell)
- Force resolution
 - 50 pN, 1 Hz – 300 kHz (mammalian)
 - 30 pN, 1 Hz – 30 kHz (turtle)
 - 10 pN, 1 Hz – 10 kHz (TRN, single molecule)
- Closed loop force control
- Operation in fluid



Overview

- PR devices
 - I-V characteristics
 - Voltage-displacement response
 - Progress towards SMFS
- PRPE devices
 - I-V characteristics
 - Voltage-force response
 - Cross-talk

The PR Device





Typical I-V Curve



Doping: 850C POCL for 15 minutes
PR: 3 μ m wide, 23 μ m long
Cantilever, 6 μ m wide, 87 μ m wide

Resistance, mobility and temperature

Thermal Breakdown




Device broke on the 3rd repeat (open circuit partway through I-V)

Runaway due to thermally generated carrier concentration becoming large.

Failure by cantilever cracking/vaporizing and aluminum vaporizing.

Noise and Power Dissipation



DC Bias, 1000x gain, 1 mHz HP filter.
4 V = 2.5 mW in PR


Thermal noise bump (returns to $1/f$ at 10 Hz). Bump frequency varies with device.

Measurements show that the noise spectrum is independent of amplifier (INA103, AD8221)


Lots of open questions. Affects the choice of optimal bias point, but not a huge problem.

Measuring PR Sensitivity





- Standard sensitivity calibration uses a piezoelectric shaker and LDV
- LDV has given me some trouble (mentioned last time)
 - Aluminum foil reduced the capacitive cross-talk 50x, should work now
- Wanted to measure the PR response more directly using AFM





PR-Only Sensitivity



Progress Towards SMFS





The PRPE Device





Piezoresistive force detection
Piezoelectric actuation

The PRPE Device



Good I-V Curves





Typical devices from a good wafer

Bad I-V Curves



Structure = 2 metal-semiconductor junctions
(electrons go from M \rightarrow S and S \rightarrow M)





Three problems:

- Exponential I-V near zero bias \rightarrow large parasitic contact resistances reduce $\Delta R/R$
- Ratio of resistances varies with voltage \rightarrow hard to balance the Wheatstone bridge
- Resistances don't match at high bias \rightarrow low Wheatstone bridge sensitivity even at high bias

Likely due to reduced surface dopant concentration from pre-AlN surface clean

Good I-V (Piezoelectric)





AlN film resistance in the gigaohm range. Higher is better.

PRPE Sensitivity


- Calibration on LDV has never worked (in contrast with the PR devices)
- From AFM testing, PR output doesn't change until the device has been visibly deformed
- After pushing some more, the PR output changes drastically (“buckles”) and behaves linearly after that

Potato Chipping (SEM)




Hypothesis: the lateral bending stiffens the cross-bar region so that force applied to the tip of the cantilever doesn't initially bend the PR region (left of the cross-bar in the first SEM)


PR+PE Buckling



PR+PE Buckling




PR+PE Buckling



PR+PE Sensitivity

Device #14 – Design #2 (400um long)




The PR output does not start to change immediately when the AFM cantilever touches the tip.

At first, all of the bending occurs out at the tip of the PR+PE device and not where the PR is located.


Eventually it reaches a threshold point, at which the device appears to controllably buckle (see images on next slide).

PR+PE Buckling


Device #13 – Design #3 (200um long)



Predicted vs. Actual




Cross-talk



Quick Cross-talk Summary

- Experimental data shown on the right
- Developed a small signal AC model. SPICE results matches within a few orders of magnitude.
- Solution: insert cross-talk signal to both sides of the Wheatstone bridge by putting all four resistors on the die



Conclusions

- PR devices
 - Work as expected
 - Working out bead attachment and functionalization in preparation for biotin-streptavidin force measurements
 - Minor issues: heating, curvature
- PRPE devices
 - Yield issues due to surface concentration reduction, about 30% usable
 - Non-linear voltage response is a show stopper, need new devices in order to be able to measure any forces
 - Cross-talk an issue for highest frequency measurements, will fix in next fab run

Next Steps

- Measure biotin-streptavidin disruption force
- Design v2 round of devices
- Fab in April

